Package 'simplex'

November 3, 2019

Title Data Reduction Software for Secondary Ion Mass Spectrometry
Version 0.1
Date 2019-11-3
Description Processes Secondary Ion Mass Spectrometery (SIMS) data within the confines of the simplex, i.e. the data space of compositions. Accommodates input files for both Cameca and SHRIMP instruments. Models the data using a combination of multinomial and logistic normal statistics. Keeps track of inter-sample error correlations caused by using a common standard for multiple samples. Includes applications for U-Pb geochronology and stable isotope geochemistry.
Author Pieter Vermeesch [aut, cre]
Maintainer Pieter Vermeesch <p.vermeesch@ucl.ac.uk></p.vermeesch@ucl.ac.uk>
Depends R (>= $3.0.0$)
Imports MASS, graphics, IsoplotR
License GPL-3
<pre>URL https://github.com/pvermees/simplex LazyData true</pre>
RoxygenNote 6.1.1
Encoding UTF-8
R topics documented:
avg_Lm
Cameca
raw_count_ratios
read_directory
subset_samples
Index

2 Cameca

avg_Lm

Average logratios

Description

Takes time resolved logratio estimates as input and produces a single spot logratio vector and covariance matrix as output.

Usage

```
avg_Lm(Lm)
```

Arguments

Lm

A time resolved list of logratio vectors and Hessian matrices.

Details

Uses a weighted mean algorithm with correlated uncertainties.

Value

A list with a single vector of logratios and its covariance matrix.

Examples

```
data(Cameca,package='simplex')
Lm <- raw_count_ratios(samp=Cameca[[1]])
aLm <- avg_Lm(Lm)</pre>
```

Cameca

Example datasets for testing simplex

Description

a Plesovice and Qinghu zircon dataset from the Cameca SIMS instrument at the Chinese Academy of Sciences Beijing.

Usage

Cameca

Format

An object of class simplex.

plot_timeresolved 3

nla	+ +	imar	eso ⁻	havi

Plot and fit raw time-resolved SIMS data.

Description

Shows the raw data for a single spot in a SIMS dataset.

Usage

```
plot_timeresolved(samp, fit = FALSE, c64 = NULL)
```

Arguments

samp	One item	of a	simplex	dataset.

fit Set to TRUE to show the maximum likelihood fit on the plot. c64 (optional) common Pb composition (206Pb/204Pb ratio).

Value

Produces a multi-panel scatter plot.

Examples

```
data(Cameca,package="simplex")
plot_timeresolved(Cameca[[1]],fit=TRUE)
```

raw_count_ratios

Time-resolved logratios

Description

Calculate the time-resolved logratios for a single spot.

Usage

```
raw_count_ratios(samp)
```

Arguments

samp

One item in an object of class simplex.

Value

A list with a vector of logratios and its covariance matrix.

read_directory

See Also

```
plot_timeresolved
```

Examples

```
data(package="simplex")
Lm <- raw_count_ratios(samp=Cameca[[1]])</pre>
```

read_directory

Read data directory

Description

Read all the input files in a data directory

Usage

```
read_directory(dname, instrument = "Cameca", suffix = NULL)
```

Arguments

dname path to the input directory

instrument text string with the type of ICP-MS. Currently only 'Cameca'.

suffix (optional) file extension of the input files.

Value

An object of class simplex, i.e. a list of lists containing the following items: ions, ions, dwelltime, detector, yield, background, cps, counts, sbm, and time.

Examples

```
datadir <- system.file(package="simplex")
dat <- read_directory(paste0(datadir,'/'),instrument='Cameca',suffix='.asc')</pre>
```

subset_samples 5

	-
SUNSAT	samples

Subset a dataset of class simplex

Description

Select a subset of samples aor standards from a simplex dataset.

Usage

```
subset_samples(dat, prefix = "Plesovice")
```

Arguments

dat A dataset to subset prefix Text string to match

Value

an object with the same class as dat

See Also

```
read_directory
```

Examples

```
data(Cameca,package="simplex")
stand <- subset_samples(Cameca,prefix='Plesovice')</pre>
```

Index